

# POPULATION ASSAY: SPORE STRIPS/DISC/THREAD/METAL

Lot # \_\_\_\_\_ POP LEVEL \_\_\_\_\_ TSA Lot # \_\_\_\_\_

CARRIER (circle one): Spore Strip Disc Thread Metal \_\_\_\_\_ Other \_\_\_\_\_

ORGANISM(S): *B. atrophaeus* *G. stearothermophilus* Other \_\_\_\_\_

## PROCEDURE:

- 1.0 Aseptically transfer 10 spore strips/discs/threads into sterile 250 ml blender cup containing 100 ml chilled processed water. [If processing metal product, aseptically transfer 1 carrier into a water blank containing 9.9 ml sterile, processed water with 0.1 ml of Tween 80 and 1ml of 3mm sterile glass beads. Vortex for 2 minutes. Insert 10 ml tube into sonicator (38.5 – 40.5 KHz, full wave industrial stack transducer) for 10 minutes. Vortex again for 2 minutes and skip to step 4.0)]
- 2.0 Blend 3-5 minutes to a homogeneous pulp of component fibers.
- 3.0 Aseptically transfer a 10 ml aliquot from the blender cup into a sterile, screw-capped 10 ml test tube. Label each tube with lot #, temperature and length of exposure.
- 4.0 Heat shock tubes in a water bath (10 minutes at 80° - 85°C for *B. atrophaeus*, 15 minutes at 95° - 100°C for *G. stearothermophilus*.) Immediately cool tubes in a water bath of 0° - 4°C.

Start Time/Temperature: \_\_\_\_\_ / \_\_\_\_\_ °C End Time: \_\_\_\_\_

Initial and Date: \_\_\_\_\_ / \_\_\_\_\_

- 5.0 Vortex the tubes for 15-20 seconds.
- 6.0 Perform serial dilutions by pipetting out 1.0 ml of the aliquot into another sterile, screwcapped 10 ml test tube containing 9.0 ml of sterile, processed water. Repeat from step 5 until desired dilution factor is reached.
- 7.0 At the dilution factors expected to yield 10-300 CFU, pipette out 1.0 ml into each of three petri plates. Repeat for final dilutions.
- 8.0 Within 20 minutes, add to each plate approximately 20 ml of TSA, pre-sterilized and cooled to 47° ± 2°C. Swirl to distribute spores evenly in agar and allow to solidify.

TSA Temperature: \_\_\_\_\_ °C Initial and Date: \_\_\_\_\_ / \_\_\_\_\_

- 9.0 Invert and incubate the plates (30° - 35°C for *B. atrophaeus*/*B. pumilus* and other mesophiles, 55° - 60°C for *G. stearothermophilus*).

Incubation Start Time/Initial & Date: \_\_\_\_\_ / \_\_\_\_\_ Incubator # \_\_\_\_\_

- 10.0 Examine all plates at 24 (±1) hours. Record on the back the number of colony forming units (CFU's) per plate. Record the average on the following page.
- 11.0 Calculate the average number of CFU's per carrier from the above data using the formulas on the following page:

Performed By: \_\_\_\_\_ Date: \_\_\_\_\_

# POPULATION ASSAY: SPORE STRIPS/DISC/THREAD/METAL

Total @ 24 hrs / number of plates counted x DF = CFU/spore carrier  
DF= Dilution factor (absolute value of the reciprocal of the dilution)  
AV= Average number of colonies per spore carrier

Incubation End Time/Initial & Date: \_\_\_\_\_ / \_\_\_\_\_

---

## CFU COUNTS AT 24 HOURS

# dilutions \_\_\_\_\_

### **24hrs**

Plates 1. \_\_\_\_\_ 2. \_\_\_\_\_ 3. \_\_\_\_\_ Total @ 24hours: \_\_\_\_\_

Total @ 24 hrs \_\_\_\_\_ / 3 x \_\_\_\_\_ (DF) = \_\_\_\_\_ (AV)CFU/Spore carrier

---

## CFU COUNTS AT 24 HOURS

# dilutions \_\_\_\_\_

### **24hrs**

Plates 1. \_\_\_\_\_ 2. \_\_\_\_\_ 3. \_\_\_\_\_ Total @ 24hours: \_\_\_\_\_

Total @ 24 hrs \_\_\_\_\_ / 3 x \_\_\_\_\_ (DF) = \_\_\_\_\_ (AV)CFU/Spore carrier

---

# of Dilutions = Dilution Factor

1 = 10

2 = 100

3 = 1000

4 = 10000

5 = 100000

6 = 1000000

Sum of the AV of both dilution / 2 =CFU/ Spore carrier

\_\_\_\_\_ / 2 =

\_\_\_\_\_ x10<sup>\_\_\_\_\_</sup> CFU/Spore Carrier

Read By: \_\_\_\_\_ Date: \_\_\_\_\_